

Anti-Gnat3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - A8Mtj3

Item Number: CM0021521



Introduction

Rabbit polyclonal antibody targeting human GNAT3, validated for WB, IF/ICC, ELISA. Recognizes gustducin alpha-3 chain involved in taste transduction and gut sugar sensing. Cross-reacts with human and mouse.

[Learn More](#)

Parameter	Value
Product Name	GNAT3 Rabbit pAb
Remarks/Alias	GDCA; HG1E; GNAT3
Species	Human
Gene ID (Human)	346562
Gene ID	346562
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-130 of human GNAT3 (NP_001095856.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	A8MTJ3
Protein Weight	40kDa
Shipping	Ice bag