

Anti-Gnaz Polyclonal Antibody For Wb, Elisa - P19086

Item Number: CM0017732



Introduction

Rabbit polyclonal antibody targeting human GNAZ (Gz-alpha), a transducer in transmembrane signaling. Validated for WB and ELISA, cross-reacts with mouse and rat. Expected band ~41 kDa.

[Learn More](#)

Parameter	Value
Product Name	GNAZ Rabbit pAb
Remarks/Alias	HG1H; gz-alpha; GNAZ
Species	Human
GeneID (Human)	2781
GeneID	2781
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human GNAZ (NP_002064.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19086
Protein Weight	41kDa
Shipping	Ice bag