


Anti-Gnb1 Polyclonal Antibody For Wb, If/Icc, Elisa - P62873

Item Number: CM0023067



Introduction

Rabbit polyclonal antibody targeting GNB1, validated for WB, IF/ICC, and ELISA. Cross-reacts with human and mouse. GNB1 is a G protein beta subunit essential for transmembrane signaling, making this antibody ideal for signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	GNB1 Rabbit pAb
Remarks/Alias	MDS; HG2A; MRD42; GNB1
Species	Human
GeneID (human)	2782
GeneID	2782
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-340 of human GNB1 (NP_002065.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P62873
Protein Weight	37kDa
Shipping	Ice bag