

Anti-Gnb1L Polyclonal Antibody For Wb, Elisa - Q9Byb4

Item Number: CM0025247

Introduction

Rabbit polyclonal antibody targeting human GNB1L, a key regulator of DNA damage response. Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for DDR signaling research.

[Learn More](#)

Parameter	Value
Product Name	GNB1L Rabbit pAb
Remarks/Alias	GY2; FKSG1; WDR14; WDVCF; DGCRK3; GNB1L
Species	Human
Gene ID (Human)	54584
Gene ID	54584
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-327 of human GNB1L (NP_443730.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYB4
Protein Weight	36kDa
Shipping	Ice bag