

Anti-Gnb2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P62879

Item Number: CM0010726



Introduction

GNB2 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat GNB2, a G protein beta subunit involved in transmembrane signaling. Ideal for signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	GNB2 Rabbit mAb
Remarks/Alias	SSS4, HG2C1, NEDHYDF, GNB2
Species	Human
Gene ID (Human)	2783
Gene ID	2783
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human GNB2 (P62879)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62879
Protein Weight	37kDa
Shipping	Ice bag