

Anti-Gnb2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62879

Item Number: CM0029304



Introduction

Rabbit polyclonal antibody against human GNB2 (G protein subunit beta-2). Validated for WB, IHC-P, IF/ICC, ELISA and cross-reacts with mouse and rat. Ideal for G protein signaling research.

[Learn More](#)

Parameter	Value
Product Name	GNB2 Rabbit pAb
Remarks/Alias	SSS4; HG2C1; NEDHYDF; GNB2
Species	Human
Gene ID (Human)	2783
Gene ID	2783
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-340 of human GNB2 (NP_005264.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62879
Protein Weight	37kDa
Shipping	Ice bag