

Anti-Gnb3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P16520

Item Number: CM0015752

Introduction

High-quality rabbit polyclonal antibody against GNB3, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat GNB3, a key G-protein beta subunit in GPCR signaling. Ideal for signal transduction and retinal disease research.

[Learn More](#)

Parameter	Value
Product Name	GNB3 Rabbit pAb
Remarks/Alias	HG2D, CSNB1H, GNB3
Species	Human
Gene ID (Human)	2784
Gene ID	2784
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-340 of human GNB3 (NP_002066.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16520
Protein Weight	37kDa
Shipping	Ice bag