

Anti-Gnb5 Polyclonal Antibody For Wb, Elisa - O14775

Item Number: CM0031020



Introduction

Rabbit polyclonal antibody targeting human GNB5 (Gbeta5), validated for WB and ELISA. Suitable for neuroscience, GPCR signaling research, and cross-reacts with mouse and rat. Optimal for RGS protein activity studies.

[Learn More](#)

Parameter	Value
Product Name	GNB5 Rabbit pAb
Remarks/Alias	GB5, HG2E, IDDCA, LADCI, gbeta5, GNB5
Species	Human
Gene ID (Human)	10681
Gene ID	10681
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human GNB5 (NP_006569.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14775
Protein Weight	44kDa
Shipping	Ice bag