


Anti-Gne Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y223

Item Number: CM0024995



Introduction

High-quality rabbit polyclonal antibody against GNE (bifunctional UDP-N-acetylglucosamine 2-epimerase/N-acetylmannosamine kinase). Validated for WB, IHC-P, ELISA, cross-reacts with human, mouse, rat. Ideal for sialic acid biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	GNE Rabbit pAb
Remarks/Alias	NM; DMRV; IBM2; Uae1; GLCNE; GNE
Species	Human
GeneID (Human)	10020
GeneID	10020
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human GNE (NP_001121699.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y223
Protein Weight	79kDa
Shipping	Ice bag