

Anti-Gng3 Polyclonal Antibody For Wb, Elisa - P63215

Item Number: CM0025696

Introduction

CamelBio offers Anti-GNG3 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human GNG3 with cross-reactivity to mouse and rat. Ideal for membrane signaling and G protein research.

[Learn More](#)

Parameter	Value
Product Name	GNG3 Rabbit pAb
Remarks/Alias	HG3D; GNG3
Species	Human
GeneID (Human)	2785
GeneID	2785
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P63215
Protein Weight	8kDa
Shipping	Ice bag