

Anti-Gng7 Rabbit Polyclonal Antibody For Wb, Elisa - O60262

Item Number: CM0014405

Introduction

High-quality rabbit polyclonal antibody against human GNG7 (HG3B), suitable for WB and ELISA. Cross-reacts with mouse and rat. Predicted 8 kDa protein. Ideal for G protein signaling research.

[Learn More](#)

Parameter	Value
Product Name	GNG7 Rabbit pAb
Alias	HG3B, GNG7
Species	Human
Gene ID	2788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-68 of human GNG7 (NP_443079.1).
Host Species	Rabbit
Category	Polyclonal Antibodies
Applications	WB, ELISA
Cross-Reactivity	Mouse, Rat
UniProt Accession	O60262
Predicted Molecular Weight	8 kDa
Shipping Condition	Ice bag