


Anti-Gngt1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P63211

Item Number: CM0018756



Introduction

Rabbit polyclonal antibody targeting GNGT1 (Guanine nucleotide-binding protein G(T) subunit gamma-T1) for WB, IF-P, ELISA. Cross-reacts with mouse and rat. Suitable for transmembrane signaling studies. P63211.

[Learn More](#)

Parameter	Value
Product Name	GNGT1 Rabbit pAb
Remarks/Alias	GNG1; HG3G1; GNGT1
Species	Human
Gene ID (Human)	2792
Gene ID	2792
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-45 of human GNGT1 (NP_068774.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P63211
Protein Weight	8kDa
Shipping	Ice bag