

Anti-Gngt2 Polyclonal Antibody For Wb And Elisa - O14610

Item Number: CM0027658

Introduction

High-quality anti-GNGT2 polyclonal antibody validated for Western blot and ELISA. Detects human and mouse GNGT2, a G protein gamma subunit involved in transmembrane signaling. Ideal for research in signal transduction.

[Learn More](#)

Parameter	Value
Product Name	GNGT2 Rabbit pAb
Remarks/Alias	GNG9; HG3I; GNGT8; G-GAMMA-8; G-GAMMA-C; GNGT2
Species	Human
GeneID (Human)	2793
GeneID	2793
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-69 of human GNGT2 (NP_113686.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O14610
Protein Weight	8kDa
Shipping	Ice bag