

Anti-Gnl1 Rabbit Polyclonal Antibody For Wb, Elisa - P36915

Item Number: CM0018945

Introduction

Rabbit polyclonal antibody against GNL1 (HSR1), raised with recombinant protein. Validated for Western blotting and ELISA; cross-reacts with human, mouse, and rat. Suitable for studying GTP-binding proteins and histocompatibility.

[Learn More](#)

Parameter	Value
Product Name	GNL1 Rabbit pAb
Remarks/Alias	HSR1, GNL1
Species	Human
GeneID (Human)	2794
GeneID	2794
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human GNL1 (NP_005266.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36915
Protein Weight	69kDa
Shipping	Ice bag