


Anti-Gnl3 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q9Bvp2

Item Number: CM0030755



Introduction

Rabbit polyclonal antibody against human GNL3 (Nucleostemin), validated for WB, IF/ICC, IP, and ELISA. Detects a 62kDa nucleolar GTP-binding protein implicated in stem cell maintenance and cancer.

[Learn More](#)

Parameter	Value
Product Name	GNL3 Rabbit pAb
Remarks/Alias	NS; E2IG3; NNP47; C77032; GNL3
Species	Human
Gene ID	26354
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human GNL3 (NP_996562.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	Q9BVP2
Protein Weight	62kDa
Shipping	Ice bag