


Anti-Gnmt Polyclonal Antibody For Wb And Elisa - Q14749

Item Number: CM0025267



Introduction

Rabbit polyclonal antibody targeting human GNMT, suitable for western blotting and ELISA. Cross-reacts with mouse and rat. Ideal for research on glycine methylation and liver metabolism.

[Learn More](#)

Parameter	Value
Product Name	GNMT Rabbit pAb
Remarks/Alias	HEL-S-182mP; GNMT
Species	Human
GeneID (Human)	27232
GeneID	27232
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-295 of human GNMT (NP_061833.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14749
Protein Weight	33kDa
Shipping	Ice bag