


## Anti-Gnpda1 Rabbit Polyclonal Antibody For Wb, Elisa - P46926

Item Number: CM0018349



### Introduction

Rabbit polyclonal antibody targeting glucosamine-6-phosphate deaminase 1, suitable for WB and ELISA applications. Detects human, mouse, and rat GNPDA1; raised against N-terminal recombinant protein. Research use for hexosamine pathway and hyaluronan synthesis studies.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | GNPDA1 Rabbit pAb                                                                                                |
| Remarks/Alias    | GPI, HLN, GNP1, GNPI, GNPDA, GNPDA1                                                                              |
| Species          | Human                                                                                                            |
| GeneID (human)   | 10007                                                                                                            |
| GeneID           | 10007                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human GNPDA1 (NP_005462.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | P46926                                                                                                           |
| Protein Weight   | 33kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |