

Anti-Gnpnat1 Rabbit Polyclonal Antibody For Wb - Q96Ek6

Item Number: CM0018837



Introduction

Rabbit polyclonal antibody against human GNPAT1 (Glucosamine 6-phosphate N-acetyltransferase), validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for Golgi and endosomal membrane studies.

[Learn More](#)

Parameter	Value
Product Name	GNPNAT1 Rabbit pAb
Remarks/Alias	GNA1; GNPAT; Gpnat1; RHZDAN; GNPAT1
Species	Human
Gene ID (Human)	64841
Gene ID	64841
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-91 of human GNPAT1 (NP_932332.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96EK6
Protein Weight	21kDa
Shipping	Ice bag