


Anti-Gnrh1 Polyclonal Antibody For Wb, If/Icc, Elisa - P01148

Item Number: CM0029995



Introduction

Rabbit polyclonal antibody targeting human GNRH1 (progonadoliberin-1), suitable for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Supports reproductive endocrinology research.

[Learn More](#)

Parameter	Value
Product Name	GNRH1 Rabbit pAb
Remarks/Alias	GRH; GNRH; LHRH; LNRH; GNRH1
Species	Human
GeneID (Human)	2796
GeneID	2796
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-92 of human GNRH1 (NP_001076580.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01148
Protein Weight	10kDa
Shipping	Ice bag