

Anti-Gnrh2 Rabbit Polyclonal Antibody For Wb, Elisa - 043555

Item Number: CM0030137

Introduction

Rabbit polyclonal antibody targeting human GNRH2 (Progonadoliberin-2), validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for reproductive biology and hormone research.

[Learn More](#)

Parameter	Value
Product Name	GNRH2 Rabbit pAb
Remarks/Alias	GnRH-II; LH-RHII; GNRH2
Species	Human
GeneID (human)	2797
GeneID	2797
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-120 of human GNRH2 (NP_001492.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43555
Protein Weight	13kDa
Shipping	Ice bag