

Anti-Gnrhr Rabbit Polyclonal Antibody For Wb, Elisa - P30968

Item Number: CM0015561



Introduction

High-quality rabbit polyclonal antibody against human GNRHR (Gonadotropin-releasing hormone receptor). Validated for Western blot (WB) and ELISA applications. Reacts with human, mouse, and rat samples. Ideal for reproductive endocrinology research.

[Learn More](#)

Parameter	Value
Product Name	GNRHR Rabbit pAb
Remarks/Alias	HH7; GRHR; LRHR; LHRHR; GNRHR1; GNRHR
Species	Human
Gene ID	2798 (Human)
Immunogen	Recombinant fusion protein containing a sequence within amino acids 1-80 of human GNRHR (NP_000397.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30968
Protein Weight	38kDa
Shipping	Ice bag