

Anti-Gnt-V/Mgat5 Rabbit Polyclonal Antibody For Wb And Elisa - Q09328

Item Number: CM0013663



Introduction

Rabbit polyclonal antibody targeting GNT-V/MGAT5 (GlcNAc-T V), validated for WB and ELISA. Detects human and mouse species. Ideal for N-glycan biosynthesis, cell signaling, adhesion, and migration research.

[Learn More](#)

Parameter	Value
Product Name	GNT-V/MGAT5 Rabbit pAb
Remarks/Alias	GNT-V; GNT-VA; MGAT5A; glcNAc-T V; GNT-V/MGAT5
Species	Human
Gene ID (Human)	4249
Gene ID	4249
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 627-741 of human GNT-V/MGAT5 (NP_002401.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q09328
Protein Weight	85kDa
Shipping	Ice bag