

Anti-Golga1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q92805

Item Number: CM0017044

Introduction

Rabbit polyclonal antibody against human GOLGA1 (golgin-97). Validated for WB, IF/ICC, ELISA; reacts with human, mouse, rat. Ideal for Golgi trafficking and E-cadherin transport studies.

[Learn More](#)

Parameter	Value
Product Name	GOLGA1 Rabbit pAb
Remarks/Alias	golgin-97; GOLGA1
Species	Human
Gene ID (Human)	2800
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human GOLGA1 (NP_002068.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92805
Protein Weight	88kDa
Shipping	Ice bag