


Anti-Golga4 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13439

Item Number: CM0007456



Introduction

Anti-GOLGA4 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets trans-Golgi network protein involved in vesicular trafficking. Ideal for Golgi apparatus research.

[Learn More](#)

Parameter	Value
Product Name	GOLGA4 Rabbit mAb
Remarks/Alias	GCP2; GOLG; p230; CRPF46; MU-RMS-40.18; GOLGA4
Species	Human
GeneID (Human)	2803
GeneID	2803
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human GOLGA4. (NP_002069.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13439
Protein Weight	261kDa
Shipping	Ice bag