


Anti-Golga4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13439

Item Number: CM0015775



Introduction

Rabbit polyclonal antibody targeting human GOLGA4 (golgin-245, p230). Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse, rat. Ideal for Golgi trafficking research.

[Learn More](#)

Parameter	Value
Product Name	GOLGA4 Rabbit pAb
Remarks/Alias	GCP2; GOLG; p230; CRPF46; MU-RMS-40.18; GOLGA4
Species	Human
Gene ID (Human)	2803
Gene ID	2803
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human GOLGA4 (NP_002069.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13439
Protein Weight	261kDa
Shipping	Ice bag