

Anti-Golga5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Tba6

Item Number: CM0017993

Introduction

Rabbit polyclonal antibody against human GOLGA5 (Golgin-84), a Golgi membrane protein involved in stack formation and retrograde transport. Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. For Golgi apparatus research.

[Learn More](#)

Parameter	Value
Product Name	GOLGA5 Rabbit pAb
Remarks/Alias	RFG5; GOLIM5; ret-II; GOLGA5
Species	Human
GeneID (Human)	9950
GeneID	9950
Immunogen	Recombinant Protein (aa 1-300)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TBA6
Protein Weight	83kDa
Shipping	Ice bag