

Anti-Golga7 Rabbit Polyclonal Antibody For Wb, Elisa - Q7Z5G4

Item Number: CM0023848



Introduction

Rabbit polyclonal antibody targeting human GOLGA7 (GCP16) for Western blot and ELISA applications. Detects both human and mouse GOLGA7, supporting research on Golgi transport and palmitoylation.

[Learn More](#)

Parameter	Value
Product Name	GOLGA7 Rabbit pAb
Remarks/Alias	GCP16; GOLGA7A; HSPC041; GOLGA3AP1; GOLGA7
Species	Human
Gene ID	51125
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-134 of human GOLGA7 (NP_001167595.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7Z5G4
Protein Weight	16kDa
Shipping	Ice bag