

Anti-Golgin-45 (Blzf1) Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H2G9

Item Number: CM0030613

Introduction

Polyclonal antibody targeting golgin-45 (BLZF1), validated for WB, IF/ICC, ELISA; reacts with human and mouse. Essential for Golgi structure and ER-to-Golgi transport, ideal for cell biology research.

[Learn More](#)

Parameter	Value
Product Name	BLZF1 Rabbit pAb
Remarks/Alias	JEM1, JEM-1, JEM-1s, GOLGIN-45, BLZF1
Species	Human
GeneID (Human)	8548
GeneID	8548
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-400 of human BLZF1 (NP_003657.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H2G9
Protein Weight	45kDa
Shipping	Ice bag