

Anti-Golph2 Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q8Nbj4

Item Number: CM0006906



Introduction

Rabbit monoclonal GOLPH2 antibody validated for WB, IF, IHC-P, ELISA. Recognizes human, mouse, rat GP73/GOLM1, a Golgi membrane protein upregulated in viral hepatitis. Suitable for research on liver disease and epithelial biology.

[Learn More](#)

Parameter	Value
Product Name	GOLPH2 Rabbit PolymAb®
Remarks/Alias	GP73; HEL46; GOLPH2; C9orf155; PSEC0257; bA379P1.3
Species	Human
Gene ID (Human)	51280
Gene ID	51280
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-350 of human GOLPH2 (NP_057632.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NBJ4
Protein Weight	45kDa
Shipping	Ice bag