

Anti-Golph2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Nbj4

Item Number: CM0013483



Introduction

Rabbit polyclonal antibody targeting human GOLPH2 (GP73, GOLM1) validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for liver disease and Golgi apparatus research.

[Learn More](#)

Parameter	Value
Product Name	GOLPH2 Rabbit pAb
Remarks/Alias	GP73; HEL46; GOLPH2; C9orf155; PSEC0257; bA379P1.3
Species	Human
GeneID (Human)	51280
GeneID	51280
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 112-401 of human GOLPH2 (NP_057632.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human,Mouse,Rat
SWISS	Q8NBJ4
Protein Weight	45kDa
Shipping	Ice bag