

Anti-Golph2 Rabbit Monoclonal Antibody For Wb, Elisa - Q8Nbj4

Item Number: CM0002929



Introduction

Rabbit monoclonal antibody against human GOLPH2 (GP73), suitable for Western blot and ELISA. Cross-reacts with mouse. Target protein weight 45 kDa, SWISS Q8NBj4. Ideal for studying Golgi function and viral liver disease biomarker research.

[Learn More](#)

Parameter	Value
Product Name	GOLPH2 Rabbit mAb
Remarks/Alias	GP73; HEL46; GOLPH2; C9orf155; PSEC0257; bA379P1.3
Species	Human
GeneID (human)	51280
GeneID	51280
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human GOLPH2 (Q8NBj4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NBj4
Protein Weight	45 kDa
Shipping	Ice bag