

Anti-Golph2/Gp73 Rabbit Polyclonal Antibody For Wb - Q8Nbj4

Item Number: CM0027488



Introduction

Rabbit polyclonal antibody against human GOLPH2/GP73, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for Golgi research and liver disease biomarker studies.

[Learn More](#)

Parameter	Value
Product Name	GOLPH2 Rabbit pAb
Remarks/Alias	GP73; HEL46; GOLPH2; C9orf155; PSEC0257; bA379P1.3
Species	Human
Gene ID (Human)	51280
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 112-401 of human GOLPH2 (NP_057632.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NBJ4
Protein Weight	45kDa
Shipping	Ice bag