


Anti-Golph3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H4A6

Item Number: CM0018637



Introduction

Rabbit polyclonal antibody targeting human GOLPH3, validated in WB, IF/ICC, ELISA; cross-reactive with mouse and rat. Suitable for studying Golgi membrane trafficking, cell migration, and mTOR signaling.

[Learn More](#)

Parameter	Specification
Product Name	GOLPH3 Rabbit pAb
Remarks/Alias	GOPP1, GPP34, MIDAS, Vps74, GOLPH3
Species	Human
GeneID (human)	64083
GeneID	64083
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-298 of human GOLPH3 (NP_071413.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H4A6
Protein Weight	34kDa
Shipping	Ice bag