

Anti-Golph4 Monoclonal Antibody For Wb, Elisa - 000461

Item Number: CM0003766

Introduction

This high-quality rabbit monoclonal antibody specifically detects GOLPH4 (GPP130) in human samples. Validated for Western blot (WB) and ELISA, it is ideal for studying endosome-to-Golgi trafficking.

[Learn More](#)

Parameter	Value
Product Name	GOLPH4 Rabbit mAb
Remarks/Alias	P138; GIMPC; GOLPH4; GPP130
Species	Human
GeneID (Human)	27333
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GOLPH4 (000461).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O00461
Protein Weight	82kDa
Shipping	Ice bag