

Anti-Gorab Polyclonal Antibody For Wb, Elisa - Q5T7V8

Item Number: CM0023053

Introduction

GORAB Rabbit Polyclonal Antibody validated for WB and ELISA. Targets human and mouse GORAB, a Golgi protein involved in COPI-mediated trafficking and glycosylation. Ideal for IVD and research applications.

[Learn More](#)

Parameter	Value
Product Name	GORAB Rabbit pAb
Remarks/Alias	GO; NTKLBP1; SCYL1BP1; GORAB
Species	Human
Gene ID	92344
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GORAB (NP_001139511.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5T7V8
Protein Weight	42kDa
Shipping	Ice bag