

Anti-Got1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P17174

Item Number: CM0027481



Introduction

Rabbit polyclonal antibody against human GOT1 (cytoplasmic aspartate aminotransferase), suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. ~46 kDa. Ideal for studying glutamate metabolism, neuroprotection, and liver function.

[Learn More](#)

Parameter	Value
Product Name	GOT1 Rabbit pAb
Remarks/Alias	AST; AST1; SGOT; cCAT; GIG18; cAspAT; ASTQTL1; GOT1
Species	Human
Gene ID (Human)	2805
Gene ID	2805
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 244-413 of human GOT1 (NP_002070.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17174
Protein Weight	46kDa
Shipping	Ice bag