

Anti-Got2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P00505

Item Number: CM0003425



Introduction

High-quality rabbit monoclonal antibody against human GOT2 (aspartate aminotransferase, mitochondrial). Validated for WB, IHC-P, and ELISA applications with cross-reactivity to mouse and rat. Ideal for metabolic and neurological research.

[Learn More](#)

Parameter	Value
Product Name	GOT2 Rabbit mAb
Remarks/Alias	KAT4; DEE82; KATIV; KYAT4; mitAAT; GOT2
Species	Human
Gene ID (Human)	2806
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GOT2 (P00505).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00505
Protein Weight	48kDa
Shipping	Ice bag