

Anti-Got2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P00505

Item Number: CM0025485



Introduction

Rabbit polyclonal antibody against human GOT2 (mitochondrial aspartate aminotransferase, kynurenine aminotransferase 4). Validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for metabolic, neurological, and fatty acid transport studies.

[Learn More](#)

Parameter	Specification
Product Name	GOT2 Rabbit pAb
Remarks/Alias	KAT4; DEE82; KATIV; KYAT4; mitAAT; GOT2
Species	Human
Gene ID (Human)	2806
Gene ID	2806
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-200 of human GOT2 (NP_002071.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00505
Protein Weight	48kDa
Shipping	Ice bag