

Anti-Gp2 Rabbit Polyclonal Antibody For Wb And Elisa - P55259

Item Number: CM0027546



Introduction

Rabbit polyclonal antibody against human GP2 (ZAP75), validated for WB and ELISA. Suitable for detecting GP2 involved in mucosal immune response as an M-cell transcytotic receptor.

[Learn More](#)

Parameter	Specification
Product Name	GP2 Rabbit pAb
Remarks/Alias	ZAP75; GP2
Species	Human
Gene ID (Human)	2813
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 111-387 of human GP2 (NP_001007243.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P55259
Protein Weight	59kDa
Shipping	Ice bag