

Anti-Gp2 Rabbit Polyclonal Antibody For Wb, Elisa - P55259

Item Number: CM0014001



Introduction

Rabbit polyclonal antibody targeting human GP2 (ZAP75), validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for studying M-cell transcytosis and mucosal immunity.

[Learn More](#)

Parameter	Value
Product Name	GP2 Rabbit pAb
Remarks/Alias	ZAP75; GP2
Species	Human
Gene ID	2813
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 111-387 of human GP2 (NP_001007243.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P55259
Protein Weight	59kDa
Shipping	Ice bag