

Anti-Gp9 Polyclonal Antibody For Wb, Ihc-P, Elisa - P14770

Item Number: CM0029156



Introduction

Rabbit polyclonal antibody against human GP9 (GPIX/CD42a), validated for WB, IHC-P, ELISA. Detects human, mouse, rat GP9, a platelet glycoprotein essential for vWF-mediated adhesion in hemostasis research.

[Learn More](#)

Parameter	Value
Product Name	GP9 Rabbit pAb
Remarks/Alias	GPIX; CD42a; GP9
Species	Human
Gene ID	2815
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-147 of human GP9 (NP_000165.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14770
Protein Weight	19kDa
Shipping	Ice bag