

Anti-Gpa33 Rabbit Monoclonal Antibody - Q99795

Item Number: CM0011626

Introduction

Rabbit monoclonal antibody targeting human GPA33 (Cell surface A33 antigen). Validated for WB, IF, IHC, ELISA. Cross-reacts with mouse and rat. Suitable for gastrointestinal and colon cancer research. Research use only.

[Learn More](#)

Parameter	Value
Product Name	GPA33 Rabbit mAb
Remarks/Alias	A33; GPA33
Species	Human
Gene ID (Human)	10223
Gene ID	10223
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99795
Protein Weight	36kDa
Shipping	Ice bag