

Anti-Gpam Polyclonal Antibody For Wb, Elisa - Q9Hcl2

Item Number: CM0019826



Introduction

Rabbit polyclonal antibody targeting human GPAM (glycerol-3-phosphate acyltransferase 1). Suitable for WB and ELISA; cross-reacts with mouse and rat. Protein weight 93kDa. Swiss-Prot Q9HCL2.

[Learn More](#)

Parameter	Value
Product Name	GPAM Rabbit pAb
Remarks/Alias	GPAM; GPAT; GPAT1; glycerol-3-phosphate acyltransferase; mitochondrial
Species	Human
Gene ID (Human)	57678
Gene ID	57678
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9HCL2
Protein Weight	93kDa
Shipping	Ice bag