


Anti-Gpc1 Polyclonal Antibody For Wb And Elisa - P35052

Item Number: CM0017773



Introduction

Rabbit polyclonal antibody against human GPC1 (Glypican-1), suitable for WB and ELISA; cross-reacts with mouse and rat. Ideal for studying heparan sulfate proteoglycan in cell signaling, prion biology, and myelination.

[Learn More](#)

Parameter	Value
Product Name	GPC1 Rabbit pAb
Remarks/Alias	glypican; GPC1
Species	Human
GeneID (Human)	2817
GeneID	2817
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-420 of human GPC1 (NP_002072.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35052
Protein Weight	62kDa
Shipping	Ice bag