

## Anti-Gpc6 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y625

Item Number: CM0028065

### Introduction

This GPC6 rabbit polyclonal antibody (pAb) is validated for WB, IF/ICC, and ELISA. It recognizes human and mouse GPC6, a heparan sulfate proteoglycan involved in WNT5A signaling, with a molecular weight of 63 kDa.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | GPC6 Rabbit pAb                                                                                                    |
| Remarks/Alias    | OMIMD1; GPC6                                                                                                       |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 10082                                                                                                              |
| GeneID           | 10082                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 320-530 of human GPC6 (NP_005699.1). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IF/ICC, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse                                                                                                       |
| SWISS            | Q9Y625                                                                                                             |
| Protein Weight   | 63kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |