

Anti-Gpd1 Rabbit Polyclonal Antibody For Wb And Elisa - P21695

Item Number: CM0027572



Introduction

GPD1 rabbit polyclonal antibody validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for studying glycerol metabolism and liver function in research applications.

[Learn More](#)

Parameter	Specification
Product Name	GPD1 Rabbit pAb
Remarks/Alias	GPD-C; HTGTI; GPDH-C; GPD1
Species	Human
GeneID (Human)	2819
GeneID	2819
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-349 of human GPD1 (NP_005267.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21695
Protein Weight	38kDa
Shipping	Ice bag