

Anti-Gpd2 Polyclonal Antibody For Wb, Elisa - P43304

Item Number: CM0026949

Introduction

Rabbit polyclonal antibody against human GPD2 (mitochondrial glycerol-3-phosphate dehydrogenase). Verified for WB & ELISA. Cross-reacts with mouse/rat. Suitable for metabolic disease and mitochondrial studies.

[Learn More](#)

Parameter	Value
Product Name	GPD2 Rabbit pAb
Remarks/Alias	GDH2; GPDH; mGDH; mGPDH; GPD2
Species	Human
GeneID (Human)	2820
GeneID	2820
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-325 of human GPD2 (NP_000399.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43304
Protein Weight	81kDa
Shipping	Ice bag