

Anti-Gpd2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P43304

Item Number: CM0009446

Introduction

CamelBio provides high-quality Anti-GPD2 Rabbit Monoclonal Antibody (P43304) validated for WB, IHC-P, ELISA. Detects human and rat GPD2, a mitochondrial glycerol-3-phosphate dehydrogenase crucial for pancreatic beta-cell glucose sensing.

[Learn More](#)

Parameter	Value
Product Name	GPD2 Rabbit mAb
Remarks/Alias	GDH2; GPDM; mGDH; mGPDH; GPD2
Species	Human
Gene ID	2820 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-325 of human GPD2 (NP_000399.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P43304
Protein Weight	81kDa
Shipping	Ice bag