

## Anti-Gpd2 Rabbit Polyclonal Antibody For Wb, Elisa - P43304

Item Number: CM0023305

### Introduction

Rabbit polyclonal antibody against human GPD2 (mitochondrial glycerol-3-phosphate dehydrogenase) for WB and ELISA applications. Suitable for human, mouse, and rat samples. Key component of pancreatic beta-cell glucose sensing.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | Anti-GPD2 Rabbit Polyclonal Antibody                                                                             |
| Remarks/Alias    | GDH2, GPDM, mGDH, mGPDH, GPD2                                                                                    |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 2820                                                                                                             |
| Gene ID          | 2820                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 25-325 of human GPD2 (NP_000399.3) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | P43304                                                                                                           |
| Protein Weight   | 81kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |