

Anti-GpIX/Cd42A/Gp9 Rabbit Monoclonal Antibody For Wb, Elisa - P14770

Item Number: CM0006427



Introduction

Rabbit monoclonal antibody targeting human platelet glycoprotein IX (GPIX/CD42a/GP9), suitable for WB and ELISA. Essential tool for hemostasis and platelet adhesion research.

[Learn More](#)

Parameter	Value
Product Name	GPIX/CD42a/GP9 Rabbit mAb
Remarks/Alias	GPIX; CD42a; GPIX/CD42a/GP9
Species	Human
GeneID	2815
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-140 of human GPIX/CD42a/GP9 (NP_000165.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P14770
Protein Weight	19kDa
Shipping	Ice bag